


Cannulated Screws / Staples

Lochschrauben / Staples

Tornillos canulados / Grapas


Viti cannulate / Staples

Vis perforée / Staples




**D-127-03 Cannulated Screws 3.5 mm and Instrument Set
196-045-17 Aluminium Case
196-127-05 Upper Tray
196-127-06 Lower Tray**

| | | | Pieces |
|-----------------------------|---|---------|--------|
| 127-035-10 To 127-035-24 | Cannulated Screw 3.5 mm (half thread) 10-24 mm | 2 each. | 16 |
| 127-035-26 To 127-035-40 | Cannulated Screw 3.5 mm (half thread) 26-40 mm | 4 each. | 32 |
| 127-035-42 To 127-035-50 | Cannulated Screw 3.5 mm (half thread) 42-50 mm | 2 each. | 10 |
| 127-045-10 To 127-045-24 | Cannulated Screw 3.5 mm (full thread) 10-24 mm | 2 each. | 16 |
| 127-045-26 To 127-045-40 | Cannulated Screw 3.5 mm (full thread) 26-40 mm | 4 each. | 32 |
| 127-045-42 To 127-045-50 | Cannulated Screw 3.5 mm (full thread) 42-50 mm | 2 each. | 10 |
| 127-085-13 | Washer 3.5 mm | | 6 |
| 127-160-27 | Cannulated Drill Bit 2.7 mm Ø | | 1 |
| 127-160-35 | Cannulated Drill Bit 3.5 mm Ø | | 1 |
| 127-168-35 | Cannulated Tap 3.5 mm Ø | | 1 |
| 127-125-10 | Guide Wire Ø 1.0 mm | | 1 |
| 127-175-35 | Cannulated screwdriver | | 1 |
| 100-072-02 | Screw Holding Forceps | | 10 |
| 127-212-35 | Tap Sleeve 2.7 mm, drill sleeve 1.0 mm | | 1 |
| 127-200-35 | Tap Sleeve 3.5 mm, drill sleeve 2.5 mm | | 1 |




D-127-04 Cannulated Screws 4.5 mm and Instrument Set
196-045-17 Aluminium Case
196-127-08 Upper Tray
196-127-09 Lower Tray

| | | | Pieces |
|-----------------------------|---|---------|--------|
| 127-055-20 To 127-055-34 | Cannulated Screw 4.5 mm (half thread) 20-34 mm | 2 each. | 16 |
| 127-055-36 To 127-055-44 | Cannulated Screw 4.5 mm (half thread) 36-44 mm | 4 each. | 32 |
| 127-055-46 To 127-055-60 | Cannulated Screw 4.5 mm (half thread) 46-60 mm | 2 each. | 16 |
| 127-065-20 To 127-065-34 | Cannulated Screw 4.5 mm (full thread) 20-34 mm | 2 each. | 16 |
| 127-065-36 To 127-065-44 | Cannulated Screw 4.5 mm (full thread) 36-44 mm | 4 each. | 20 |
| 127-065-46 To 127-065-60 | Cannulated Screw 4.5 mm (full thread) 46-60 mm | 2 each. | 18 |
| 127-085-16 | Washer 4.5 mm | | 6 |
| 127-160-32 | Cannulated Drill Bit 3.2 mm \varnothing | | 1 |
| 127-160-45 | Cannulated Drill Bit 4.5 mm \varnothing | | 1 |
| 127-168-45 | Cannulated Tap 4.5 mm \varnothing | | 1 |
| 127-125-16 | Guide Wire \varnothing 1.6 mm | | 1 |
| 127-175-45 | Cannulated screwdriver | | 10 |
| 100-072-02 | Screw Holding Forceps | | 1 |
| 100-115-45 | Tap Sleeve for 4.5 mm Screw | | 1 |



D-127-07 Cannulated Screws 7.0 mm and Instrument Set
196-045-17 Aluminium Case
196-127-01 Upper Tray
196-127-02 Lower Tray

| | | | Pieces |
|-----------------------------|---|---------|--------|
| 127-070-30 To 127-070-65 | Cannulated Screw 7.0 mm \varnothing , 16 mm 30-35-40-45-50-55-60-65 mm | 2 each. | 16 |
| 127-070-70 To 127-070-92 | Cannulated Screw 7.0 mm \varnothing , 16 mm 70-75-80-85-90-95-100 mm | 4 each. | 28 |
| 127-070-93 To 127-070-96 | Cannulated Screw 7.0 mm \varnothing , 16 mm 105-110-115-120 mm | 2 each. | 8 |
| 127-080-45 To 127-080-65 | Cannulated Screw 7.0 mm \varnothing , 32 mm 45-50-55-60-65 mm | 1 each. | 5 |
| 127-080-70 To 127-080-92 | Cannulated Screw 7.0 mm \varnothing , 32 mm 70-75-80-85-90-95-100 mm | 2 each. | 14 |
| 127-080-93 To 127-080-96 | Cannulated Screw 7.0 mm \varnothing , 32 mm 105-110-115-120 mm | 1 each. | 4 |
| 127-085-19 | Washer 7.0 mm \varnothing | | 6 |
| 127-160-50 | Drill Bit 5.0 mm \varnothing | | 1 |
| 127-168-70 | Tap 7.0 mm \varnothing | | 1 |
| 127-175-70 | Screw Driver | | 1 |
| 127-227-70 | Tissue Protector Sleeve | | 1 |
| 127-125-20 | Guide Wire 2.0 mm \varnothing 230 mm | | 10 |
| 100-072-20 | Screw Holding Forceps | | 1 |



3.5 mm \varnothing with over flats S.W. 2.5 mm half / full threaded drill bit 2.7 mm

\varnothing 3.5mm, 4.5mm, 7.0mm



4.5 mm \varnothing with over flats S.W. 3.5 mm half / full threaded drill bit 3.2 mm

| Dia of thread | 3.5 mm | 3.5 mm |
|------------------|------------|------------|
| Width over flats | 2.5 mm | 2.5 mm |
| Drill Bit | 2.7 mm | 2.7 mm |
| Length | | |
| 10 mm | 127-035-10 | 127-045-20 |
| 12 mm | 127-035-12 | 127-045-22 |
| 14 mm | 127-035-14 | 127-045-24 |
| 16 mm | 127-035-16 | 127-045-26 |
| 18 mm | 127-035-18 | 127-045-28 |
| 20 mm | 127-035-20 | 127-045-30 |
| 22 mm | 127-035-22 | 127-045-32 |
| 24 mm | 127-035-24 | 127-045-34 |
| 26 mm | 127-035-26 | 127-045-36 |
| 28 mm | 127-035-28 | 127-045-38 |
| 30 mm | 127-035-30 | 127-045-40 |
| 32 mm | 127-035-32 | 127-045-42 |
| 34 mm | 127-035-34 | 127-045-44 |
| 36 mm | 127-035-36 | 127-045-46 |
| 38 mm | 127-035-38 | 127-045-48 |
| 40 mm | 127-035-40 | 127-045-50 |
| 42 mm | 127-035-42 | 127-045-52 |
| 44 mm | 127-035-44 | 127-045-54 |
| 46 mm | 127-035-46 | 127-045-56 |
| 48 mm | 127-035-48 | 127-045-58 |
| 50 mm | 127-035-50 | 127-045-60 |

| Dia of thread | 4.5 mm | 4.5 mm |
|------------------|------------|------------|
| Width over flats | 3.5 mm | 3.5 mm |
| Drill Bit | 3.2 mm | 3.2 mm |
| Length | | |
| 20 mm | 127-055-20 | 127-065-20 |
| 22 mm | 127-055-12 | 127-065-22 |
| 24 mm | 127-055-14 | 127-065-24 |
| 26 mm | 127-055-16 | 127-065-26 |
| 28 mm | 127-055-18 | 127-065-28 |
| 30 mm | 127-055-20 | 127-065-30 |
| 32 mm | 127-055-22 | 127-065-32 |
| 34 mm | 127-055-24 | 127-065-34 |
| 36 mm | 127-055-26 | 127-065-36 |
| 38 mm | 127-055-28 | 127-065-38 |
| 40 mm | 127-055-30 | 127-065-40 |
| 42 mm | 127-055-32 | 127-065-42 |
| 44 mm | 127-055-34 | 127-065-44 |
| 46 mm | 127-055-36 | 127-065-46 |
| 48 mm | 127-055-38 | 127-065-48 |
| 50 mm | 127-055-40 | 127-065-50 |
| 52 mm | 127-055-42 | 127-065-52 |
| 54 mm | 127-055-44 | 127-065-54 |
| 56 mm | 127-055-46 | 127-065-56 |
| 58 mm | 127-055-48 | 127-065-58 |
| 60 mm | 127-055-50 | 127-065-60 |



**7 mm \varnothing with over flats SW 5.0 mm
16 mm length of thread, drill bit 5.0 mm**



**7 mm \varnothing with over flats SW 5.0 mm
32 mm length of thread, drill bit 5.0 mm**

| | | | |
|------------|-------|------------|--------|
| 127-070-25 | 25 mm | 127-070-75 | 75 mm |
| 127-070-30 | 30 mm | 127-070-80 | 80 mm |
| 127-070-35 | 35 mm | 127-070-85 | 85 mm |
| 127-070-40 | 40 mm | 127-070-90 | 90 mm |
| 127-070-45 | 45 mm | 127-070-91 | 95 mm |
| 127-070-50 | 50 mm | 127-070-92 | 100 mm |
| 127-070-55 | 55 mm | 127-070-93 | 105 mm |
| 127-070-60 | 60 mm | 127-070-94 | 110 mm |
| 127-070-65 | 65 mm | 127-070-95 | 115 mm |
| 127-070-70 | 70 mm | 127-070-96 | 120 mm |

| | | | |
|------------|-------|------------|--------|
| 127-080-40 | 40 mm | 127-080-90 | 90 mm |
| 127-080-45 | 45 mm | 127-080-91 | 95 mm |
| 127-080-50 | 50 mm | 127-080-92 | 100 mm |
| 127-080-55 | 55 mm | 127-080-93 | 105 mm |
| 127-080-60 | 60 mm | 127-080-94 | 110 mm |
| 127-080-65 | 65 mm | 127-080-95 | 115 mm |
| 127-080-70 | 70 mm | 127-080-96 | 120 mm |
| 127-080-75 | 75 mm | | |
| 127-080-80 | 80 mm | | |
| 127-080-85 | 85 mm | | |




| | | |
|------------|---------------------|--------------------------------------|
| 127-085-13 | 13 mm \varnothing | Washers for Cannulated Screws 3.5 mm |
| 127-085-16 | 16 mm \varnothing | Washers for Cannulated Screws 4.5 mm |
| 127-085-19 | 19 mm \varnothing | Washers for Cannulated Screws 7.0 mm |




- 127-125-10 Guide wire 1.0 mm for Cannulated Screws 3.5 mm \varnothing 230 mm
- 127-125-16 Guide wire 1.6 mm for Cannulated Screws 4.5 mm \varnothing 230 mm
- 127-125-20 Guide wire 2.0 mm for Cannulated Screws 7.0 mm \varnothing 230 mm




- 127-160-27 Cannulated drill bit 2.7 mm for 3.5 mm Screws
- 127-160-35 Cannulated drill bit 3.5 mm for 3.5 mm Screws
- 127-160-32 Cannulated drill bit 3.2 mm for 4.5 mm Screws
- 127-160-45 Cannulated drill bit 4.5 mm for 4.5 mm Screws
- 127-160-50 Cannulated drill bit 5.0 mm for 7.0 mm Screws



- 127-168-35 Cannulated Tap for 3.5 mm Screws
- 127-168-45 Cannulated Tap for 4.5 mm Screws
- 127-168-70 Cannulated Tap for 7.0 mm Screws



- 127-173-35 Cannulated Screw Driver for 3.5 mm screws (fiber handle)
- 127-173-45 Cannulated Screw Driver for 4.5 mm screws (fiber handle)
- 127-173-70 Cannulated Screw Driver for 7.0 mm screws (fiber handle)



- 127-175-35 Cannulated Screw driver for 3.5 mm Screws with (T-handle)
- 127-175-45 Cannulated Screw driver for 4.5 mm Screws with (T-handle)
- 127-175-70 Cannulated Screw driver for 7.0 mm Screws with (T-handle)



127-200-35 Tap Sleeve 3.5 mm for 3.5 mm Screws
Drill Sleeve 2.5 mm



127-212-35 Tap Sleeve 2.7 mm for 3.5 mm Screws
Drill Sleeve 1.0 mm



100-115-08 Tap Sleeve 4.5 mm for 4.5 mm Screws
Drill Sleeve 3.2 mm



127-227-70 Tissue Protector Sleeve
for 7.0 mm Screws



127-245-20 Parallel Target Device for Guide Wires \varnothing 2.0 mm



Stainless Steel



- 128-002-16** Barbed Staples 8.7 x 16 x 2 mm
- 128-002-19** Barbed Staples 8.7 x 19 x 2 mm
- 128-002-22** Barbed Staples 8.7 x 22 x 2 mm
- 128-002-25** Barbed Staples 8.7 x 25 x 2 mm




| | Diameter | Width | Length |
|-------------------|----------|-------|--------|
| 128-015-22 | 1.5 mm | 22 mm | 19 mm |
| 128-015-16 | 1.5 mm | 16 mm | 19 mm |
| 128-025-22 | 2.5 mm | 22 mm | 19 mm |
| 128-025-16 | 2.5 mm | 16 mm | 19 mm |



128-033-48
24 x 33 x 4.8 mm

128-033-95
24 x 33 x 9.5 mm

128-033-14
24 x 33 x 14 mm



128-130-25 Staples Extractor




128-132-16 Staples Driver




128-134-19 Staples Inserter




- 128-250-06 6 mm
- 128-250-08 8 mm
- 128-250-11 11 mm
- 128-250-14 14 mm




128-270-18 Impactor




- 128-322-09 9.5 x 22 mm
- 128-322-16 16 x 22 mm
- 128-322-22 22 x 22 mm
- 128-322-28 28 x 22 mm




128-271-02 Impactor for 132-322-09 To 132-322-28
 128-271-04 Impactor for 132-340-16 To 132-355-15



| Cat. No.: | Width | Length |
|------------|-------|--------|
| 128-340-16 | 16 mm | 22 mm |
| 128-340-22 | 22 mm | 22 mm |
| 128-340-28 | 28 mm | 22 mm |



| Cat. No.: | Width A | Length B | Height C |
|------------|---------|----------|----------|
| 128-348-05 | 24 mm | 30 mm | 5 mm |
| 128-348-10 | 24 mm | 30 mm | 10 mm |
| 128-348-15 | 24 mm | 30 mm | 15 mm |



| Cat. No.: | Width | Length | Height |
|------------|-------|--------|--------|
| 128-355-05 | 24 mm | 30 mm | 5 mm |
| 128-355-10 | 24 mm | 40 mm | 10 mm |
| 128-355-15 | 24 mm | 40 mm | 15 mm |